

Clear5

PRECISION ALIGNERS

**Perfect Smile,
Perfect Precision.**

FORESTADENT® UK
GERMAN PRECISION IN ORTHODONTICS






Engineered for clinical excellence and consistent results.

Developed with a state-of-the-art multilayered material, Clear5 delivers optimal force, sustained pressure, and exceptional comfort—all while minimising the need for refinements.

From complex rotations to precise spacing, Clear5 is designed to treat a wide range of cases with accuracy and control.

Manufactured by a leading German aligner company and treatment planning designed by Orthodontists, Clear5 combines reliability, speed, and clear competitive pricing to deliver smarter, more effective aligner treatment.



Trusted by Orthodontists. Loved by Patients.



**Proven with over
500,000 Smiles**

since 2016



15.8% Refinement Rate

significantly below the industry average



7 Day Manufacture

once treatment approved



3 Million

aligners forecast to be produced this year



**Advanced Treatment
Planning Solutions**

Clear5 combines intelligent planning tools with expert insight to deliver precise outcomes. Whether you manage cases independently or collaborate with our team, you'll benefit from a system built to optimise accuracy, efficiency, and confidence at every stage.



**Transform Smiles
in as little as 5 Months**

with advanced material technology and long-lasting force retention, Clear5 aligners drive efficient, controlled tooth movement—delivering noticeable results in as little as 5 months.


Triple-Layer Structure for Precision and Performance

Clear5 aligners are constructed using a unique tri-layer material system — engineered to deliver optimal biomechanics, superior comfort, and reliable aesthetics throughout the treatment journey.



Maintains Force For Over 7 Days

Ensuring maximum plastic efficiency, guiding teeth as planned.



Lower Forces

On average, it exerts 40% lower insertion forces than conventional aligners, significantly reducing discomfort.



The Material Offers 3D Elastic Properties

Ensures precise rebound for all forces, including rotation.




8-10 Bars Pressure to Thermoform

We 3D print our models at 70 micron layer.

Layer Structure

Ensures each aligner performs with accurate force application, enhanced wear resistance and patient comfort.



Outer Shell Layer

Protective barrier

- Exceptional wear resistance against masticatory forces
- Stain-resistant for lasting transparency
- Maintains structural integrity over time

This high-durability exterior safeguards both aesthetics and performance.

Core Elastomer Layer

Force distribution

- Optimized elasticity for biomechanically efficient force delivery
- Increased flexibility for patient comfort
- Smooth performance during complex tooth movements

Designed for comfort without compromising clinical effectiveness.

Inner Precision Layer

Targeted fit and control

- High rigidity ensures tight, adaptive fit
- Enhances control over tooth movement
- Supports predictable, staged progression

Engineered for controlled and precise results.

Treatable Malocclusions

Clear5 aligners are engineered to address a wide range of malocclusion types — from minor aesthetic concerns to complex orthodontic cases. Leveraging advanced digital tools, Clear5 enables the creation of precise, personalized treatment plans that deliver reliable, controlled outcomes over time.

Crossbite



Open Bite



Crowding



Deep Bite



Un-Aligned



**Crown & Bridge
or Existing Implant**




Extractions



Displaced



Class II or III



Spacing



Excessive Overjet





Real Results

From deep bites and crossbites to Class II corrections and arch refinement, Clear5 delivers predictable outcomes with a discreet, patient-friendly experience.

Why Professionals Trust Clear5

- ✓ Treat simple to complex cases
- ✓ Aligners fabricated from advanced 3-layer aligner technology for precise, gentle movement
- ✓ Designed for comfort, clarity, and confidence throughout treatment

Workflow

Clear5 redefines the orthodontic workflow, streamlining every step from digital scans to precise 3D planning and custom aligner fabrication. The system is designed for clinical accuracy, speed, and simplicity — helping practices deliver exceptional results more efficiently.

By adopting Clear5, orthodontists can ensure high-quality, predictable outcomes while providing an improved treatment experience for patients and clinicians alike.



Intraoral Scan

Fast, Accurate Digital Impressions.

Accurate digital impressions are essential for effective aligner treatment. It's essential to generate precise 3D images of the dentition and gingiva to ensure optimal data for treatment planning. Clear5 is compatible with any intra oral scanner.

- Easy and accurate digital impression workflows
- Advanced scanning technology
- Rapid and precise data acquisition



Planning

Smarter Planning with Every Case.

Clear5's Motion Viewer 2.0 software translates scan data into comprehensive digital treatment plans. Dynamic 3D simulations models illustrate predicted tooth movements, helping clinicians visualize and refine treatment strategy with confidence.

- Full-featured orthodontic planning capabilities
- Sophisticated staging tools
- Integrated support for CBCT (DVT) data



Treatment

Efficient, Predictable, Precise.

Once your treatment plan is approved, you'll receive the Clear5 delivery kit, which contains all essential components for a smooth start. The smart packaging is well presented and includes all relevant usage and application guides for the doctor and patient.

- **Comprehensive treatment documentation**
- **Clear, step-by-step clinical instructions**
- **Streamlined, modern packaging for enhanced efficiency**



Monitoring

Real-Time Oversight.

Maintain full visibility and control of patient progress with the Clear5 Portal. This dedicated digital interface supports case submission, real-time tracking, treatment status updates, and direct messaging with the Clear5 support team — all tailored to enable optimal treatment outcomes.

- **Complete case management in one platform**
- **Live status tracking and instant updates**
- **DoctorPro online scheduling integration for streamlined workflow**



Retention

Protect Results. Prevent Relapse.


Retention is critical to maintaining treatment outcomes following active orthodontic correction. Clear5 includes a set of custom-fitted, transparent retainers engineered to preserve final tooth positions and prevent relapse.

- **Precision-fitting retainers for maximum stability**
- **Durable, long-lasting thermoplastic materials**
- **Low-maintenance and hygienic design**

Our Carbon Footprint

Forestadent is committed to the pursuit of global and environmental sustainability.

Clear5 meets these environmental goals. Our innovative poly pouches are 100% recyclable, offering a smarter, eco-conscious alternative without compromising quality or durability. It's just one of the many choices we make to reduce waste, lower emissions and put the environment first.



Treatment Packages

 <p>MICRO 0-5 months Up to 10 aligner sets Ideal for: Relapse cases post-orthodontic treatment</p> <ul style="list-style-type: none"> • Lower central incisor mild rotation relapse • Slight flaring of upper lateral incisors • Pre-veneer minor spacing adjustment • Retainer replacement with minor correction 	 <p>MINI 5-10 months Up to 20 aligner sets Ideal for: Very Mild Class I, II, III cases</p> <ul style="list-style-type: none"> • Mild relapse of lower incisors • Diastema closure (1-2 mm gap) • Minor midline shift • Minor crowding 	 <p>MAX 10 - 15 months Up to 30 aligner sets Ideal for: Mild to moderate Class I, II or III cases</p> <ul style="list-style-type: none"> • Mild lower anterior crowding • Small anterior open bite • Relapse after braces (crowding) • Class I with 4 mm overjet • Mild deep bite correction 	 <p>UNLIMITED 15+ months Unlimited aligners Ideal for: Complex Class I, II, or III cases</p> <ul style="list-style-type: none"> • Moderate to severe crowding (>4-5 mm) • Significant deep bite or open bite • Multiple posterior crossbites • Complex rotations (e.g., canines or premolars rotated >45°)
--	--	--	---

Product Options

	MICRO	MINI	MAX	UNLIMITED
Stages	10	20	30	Unlimited
Free Refinements	1	2	2	Unlimited
Expiration (Years)	1	3	3	3
Retainer Sets	1	1	1	1



forestadentuk.com

Email: info@forestadent.co.uk **Call:** 01908 227851

Forestadent Ltd, 5 Diamond Court, Opal Drive,
Fox Milne, Milton Keynes, MK15 0DU

